

HPLC Analysis of Fat-Soluble Vitamins

Application Note

AN0061

INTRODUCTION

Vitamins are critical compounds which are essential for normal metabolism function. They are naturally found in many foods but are also often added to processed food products. Additionally, vitamin supplements are a growing trend among people whose diet is restricted. Vitamins are categorised into two groups; water soluble and fat soluble. The most common fat-soluble vitamin supplements are Vitamin A (retinol), Vitamin A Acetate (Retinol Acetate), Vitamin D₂ (Ergocalciferol), Vitamin D₃ (Cholecalciferol), Vitamin E (dl- α -tocopherol), Vitamin E Acetate (dl- α -tocopherol acetate) and Vitamin K₁ (Phylloquinone).

SCION Instruments developed a qualitative method for the simultaneous identification of seven fat-soluble vitamins. As individual vitamins are unstable, it is recommended for quantitative analysis that each vitamin component is prepared and analysed with individual analytical methods.

EXPERIMENTAL

A SCION 6000 HPLC with DAD was used with a C18 reverse phase column for the simultaneous identification of seven target compounds. Utilising the capability of the DAD to select multiple wavelengths it was possible to identify all seven target analytes over three different wavelengths.

Analytical standards were prepared with a range from 0.01mg/L to 10mg/L for Vitamin A, 0.1mg/L to 100mg/L for Vitamin A acetate, Vitamin D₂, Vitamin D₃, and Vitamin K₁ whilst the calibration range of Vitamin E and Vitamin E acetate was 1mg/L to 1000mg/L. Standard dilutions were made from stocks using methanol. Samples included a Vitamin E enriched supplement and medicated eye drops. All samples were prepared with a 1:10 methanol dilution before being passed through a 0.2 μ m filter. Analytical conditions for the HPLC-DAD can be found in Table 1.

Table 1. Analytical conditions of the HPLC-DAD system

Conditions	
Column	C18 5 μ m x 4.6mm ID x 250mm
Column Temp	40°C
Mobile Phase	Acetonitrile: Methanol (60:40 v/v)
Flow Rate	1mL/min
Injection Vol	10 μ L
DAD	265nm, 280nm and 325nm

RESULTS

Excellent linearity was observed throughout the analysis with each individual calibration curve exhibiting an R² of 0.9999 or greater. Figures 1 and 2 highlight the calibration curves for Vitamin A and Vitamin E.




Figure 1. Calibration curve of Vitamin A




Figure 2. Calibration curve of Vitamin E

As shown in the above figures, it is possible to analyse a wide range of concentrations on the SCION HPLC-DAD without compromising on linearity. This eliminates the need for difficult sample preparation or method adjustments.

Table 2 details the peak identifiers of the seven fat-soluble vitamins analysed and the associated peaks observed in the chromatogram (Figure 3). Additionally, the varying wavelengths used for detection is also listed.

Table 2. Vitamin name, detector wavelength and peak number

Peak No.	Vitamin Name	Wavelength (nm)
1	Vitamin D ₂	265nm
2	Vitamin D ₃	265nm
3	Vitamin K ₁	265nm
4	Vitamin E	280nm
5	Vitamin E Acetate	280nm
6	Vitamin A	325nm
7	Vitamin A Acetate	325nm




Figure 3. Chromatogram of analytical standard

Both the eye drop sample and supplement sample were analysed with peak identification further confirmed via absorbance spectrum comparisons of both sample and standard. The ability for spectrum comparisons is a vital part of the Compass CDS software offering confidence in results.

Figures 4-8 detail the chromatograms of both the eye drop and supplement samples as well as example absorbance spectrum comparison for Vitamin A, Vitamin E Acetate and Vitamin E.




Figure 4. Chromatogram of eye drop sample




Figure 5. Absorbance spectrum comparison; Vitamin A (eye drop)




Figure 6. Absorbance spectrum comparison; Vitamin E Acetate (eye drop)




Figure 7. Chromatogram of supplement sample




Figure 8. Absorbance spectrum comparison; Vitamin E (supplement sample)

As shown by the above figures, Vitamin A and Vitamin E acetate were detected in the eye drop sample whereas only Vitamin E was detected in the supplement sample. Each identified sample peak gave the same absorbance spectrum as that in the analytical standard, giving confidence in the confirmation of analytes.

CONCLUSION

Simultaneous identification of seven fat-soluble vitamins was easily achieved using the SCION HPLC with Diode Array Detection and a reverse phase C18 column. A single method utilising multiple wavelengths allowed simultaneous identification of a wide concentration range of vitamins. Compass CDS software allows easy comparison of absorbance spectra for confirmation in identification.