APPLICATION NOTE AN136





Determination of Terpenes in Cannabis by Headspace GC-SQMS

INTRODUCTION

Hemp and marijuana are becoming more popular because of legalisation in multiple countries. The increase of this market also increases the offer of products and therefore the need to monitor different components. One of these components are cannabinoids or just short potency analysis, which can be found in our portfolio of Cannabis applications (AN091,92,93,130). This application focusses on terpenes analysis using GC-SQMS, for Terpene analysis using GC-FID see application note 133.

Terpenes are organic compounds in marijuana that provide the distinguished aroma and flavour. Every species or strain of marijuana has it's own distinguished profile and therefore it is a perfect method to analyse specific strains. These terpenes do not only provide flavour or aroma but also support other components on producing the desired effects, terpenes can help in producing a calming affect but can also give you more energy (energizing effect).

Figure 1 shows the SCION Instruments 436 GC with single quad (SQ) mass spectrometer together with the Scion HT3 Headspace sampler. On which this method is applicable to.



Figure 1. SCION HT3 Headspace Sampler together with the SCION Instruments 436 GC with SQ-MS.



Determination of Terpenes in Cannabis by GC-FID

EXPERIMENTAL

This analysis can be implemented on the 436-GC and the 456-GC platform. The analysis was performed on the Scion 456-GC analyser equipped with a single quad mass spectrometer and Scion HT3 Headspace autosampler.

There are more than one hundred terpenes known. For this application the focus is on the 21 most important components, which are shown in table 1.

Table 1: Terpenes components

Nr.	Component	Nr.	Component
1	α-Pinene	12	Terpinolene
2	Camphene	13	Linalool
3	β-Mycrene	14	Isopulegol
4	β-Pinene	15	Geraniol
5	∆-3-Carene	16	β-Caryophllene
6	α -Terpinene	17	α-Humulene
7	d-Limonene	18	Nerolidol-1
8	P-Cymene	19	Guaiol
9	Eucalyptol	20	Caryphyllene oxide
10	Ocimene	21	α-Bisabolol
11	V -Terpinene		

The focus of the calibration lies on these 21 components, the calibration curves for the terpenes standards were prepared between 0.5 and $16 \mu g/ml$.

The Quality Control (QC) sample was made from $4 \mu g/ml$ for all components.

To analyse terpenes in marijuana an extraction has to be performed. 0.1g of the sample was added to 30mL of methanol before shaking vigorously for one minute. The extract was left to settle for 30 minutes before being filtered. 2 mL was then transferred into a 20 ml Scion headspace vial, ready for injection.

Please keep in mind that the sample preparation is different than the potency preparation where the sample is heated. Table 2 details the GC, MS and headspace parameters used throughout this method.

Table 2. Instrumentation operating conditions.

Injector	Splitless 50:1, 230 °C				
Column	624-MS				
Oven Program	50°C (1.0 min), 10°C/min to 250°C (2.0 min) Helium, 1.6 ml/min				
Carrier					
Autosampler	Scion HT3				
Software	MS work station/ HT3 Teklink				
Oven temperature	150°C				
Transfer line temperature	150°C				
Sample temperature	120°C				
Sample equilibrium time	30 min				
MS transfer line temperature	250°C				
lon source temperature	200°C				
Ionization mode	EI 2.00				
Scan start					
Scan mode	Full Scan				

RESULTS

The precision of the method was obtained by seven consecutive injections containing 4 μ g/ml of each component. While reviewing the data it is important to keep in mind that these values are on trace level, and that the headspace injection can cause an extra variance in the data.

It was shown that the RSD for the components lied between 2.0 and 7.0 %., which is really good for this method. The repeatability of each component can be found in table 3.

All the chromatographic peaks were identified with the NIST library and integrated. β -Mycrene and β -Pinene have the same quantification ion and almost have the same retention times, it is advised to check the integration of these components after analysis.

Figure 2 shows examples of these target peaks from a 4 μ g/ml standard containing a mix of all the components.



α-Pinene Camphene Sample Spectrum for Scan: 2057 R BP 93 (3.239e+7=100%) rep 1 4 .xm 93 3.239e+7 93 93 93.183e+6 ep 1 4 1009 1009 4(759 75% 50% 50% 67 24 25% 25% L.Î. haller 0% 0% 300 400 m/ 200 300 400_{m/2} 200 30 100 30 100 m for .alpha.-Pinen n for Camphene an: 1 RT 93 8.902 min an: 1 RT: 93 100% 1009 75% 75% 8 20 20 8 20 20 39 385 50% 50% 25% 25% ավե 0% 300 400 m/j ectrum for Scan: 1974 RT: 8.955 ≥100%) rep 1 4 .xms 0% 200 300 400 m/z Raw Sample Spectrum for Scan. 2057 RT: 9.313 BP 93 (3.244e+7=100%) rep 14 .xms 93 3.24fe+7 100 10 10 BP 93 (9.283e+6 100% 100% 9.28 75% 75% 50% 50%-67 24 q 25% 25% 9.014=-6 hullu 0% 0% 8.9 100 200 300 400 9¹7 minutee 200 300 400 8.7 9¹3 minute 9.3 9.5 100 9.1 9.1 **β-Mycrene β-Pinene** Sample Spectrum for Scan: 2186 RT: 9.865 mi 3P 93 (1.309e+7=100%) rep 1 4 .xms Quan; 93.0 ((+) 41.0:450.0) 2177 RT: 9.827 mir Quan; 93.0 ((+) 41.0:450.0) P 93 (5.201e+7=100%) rep 1 4 .xms 1009 93 5.201e+7 1.30 e+7 50 75% 50-75% 50% 50% 25% 25% 1.285e+7 40 лH 40 0%-0%-400 m/z 300 300 400 200 100 100 200 m for .beta.-Myrce Dine 1:1 RT an: 1 RT: 9.814 r <u>aı.</u> 41 72? ∮ 1009 1009 30 30-75% 75% MCps MCps 41 400 50% 50% 25% 259 20 20 130 90 1.11 0% 0% 200 300 400 m/z Raw Sample Spectrum for Scan 2186 RT 9.865 93 93 93 93 93 93 93 93 93 93 93 93 93 93 93 1.3 the+7 300 400 100 200 _____ 300 400 m/2 Spectrum for Scan: 2177 RT: 9.827 +7=100%) rep 1 4 .xms raw Sample Sp BP 93 (5.202e+ 10 10 100% 1009 5.202e+7 75% 75% 50% 50% 25% 1.292e+7 25% 0%-0% 400 m/; 300 9.8 10.0 200 300 400 m/ 100 200 9.6 10.2 minute 100 Δ-3-Carene α -Terpinene Sample Spectrum for Scan: 2286 RT: 10.293 m BP 93 (4.107e+7=100%) rep 1 4 .xms 93 4.107e+7 93.0 ((+) 41.0:450.0) m for Scan: 2332 RT: 10.491 mi Quan; 121.0 ((+) 41.0:450.0) BP 93 (6.007e+6=100%) rep 1 4 .xm 1009 1009 6.00 30 +6 75% 75% 50% 50% 25% 9.9 25% 25 <u>multi</u> 0%-0% 400 m/ 300 400 m/z 300 200 200 100 100 um for 3-Carene for Alpha-Terpin 20 30 1 RT: 10.285 m 93 049 1 R 100% 100 75% 75% ACps MCps 15 50% 79 320 136 151 50% 20 25% 25% հեր 0% 10 0% Raw Sample Spectrum for Scan: 2286 RT: 10.293 BP 93 (4.110e+7=100%) rep 14 .xms 93 4.110e+7 200 300 400 m/: RT: 10.491 rep 1 4 .xms 3P 93 (6. 100% 1009 6.0 75% 75% 50%

9.959e+6 25%

100

200

300

400 m/

0%

10.7

10.1

10.3

10.5

9.9

Figure 2. Target compound example peaks.







50%

25%

0%

10.7

10.9

11.1

11.3

121

200

300

400

11.2

11.4

11.6

100





11.8

12.0

LILL.

100

200

300

400

0%



50%

25%

0%

18.0

18.4

18.2

18.6

18.8 minutes

ԱԱԽ

100

400 m/z

300

200



Table 3. Target compound example peaks part 3.







Table 3. Target compound example peaks part 4.

The calibration curves for the Terpene standards were prepared between 0.5 to 16 μ g/ml. Most components had a correlation coefficient (R2) of 0.999 or greater which is perfect for this application. Figure 4 shows an example of the calibration curve for three components.

20.8 minut 400

20.4

20.2

20.6

From the linearity the limit of detection (LOD) and limit of quantitation (LOQ) was calculated.

The QC sample was analysed using seven injections, the average concentrations are shown in table 3 together with the LOD and LOQ.



Caryphyllene oxide



Figure 4: Three example calibration curves of the Terpenes standard.

After analysis of the marijuana extract it showed that the sample contained α -Pinene, Camphene, β -Mycrene, β -Pinene, Linalool, α -Humulene, Guaiol and α -Bisabolol. The concentration of these components are also shown in table 3

In this app note a sample volume of 2 ml is advised, when the components are below the LOD and LOQ it is possible to increase the sample volume as long if the total volume is 10 ml or another option is enrichment of the sample.

Table 5. Results of affecting reperies analysis	Table 3.	Results of	different ⁻	Terpenes	analysis.
---	----------	-------------------	------------------------	----------	-----------

Nr.	Component	R ²	LOD µg/ml	LOQ µg/ml	QC µg/ml	Repeatability %	Sample µg/ml
1	α-Pinene	0.9996	0.24	0.71	4.1	4.	15.1
2	Camphene	0.9999	0.12	0.35	4.3	3.93	0.39
3	β-Mycrene	0.9997	0.21	0.65	4.0	3.13	2.80
4	β-Pinene	0.9997	0.21	0.65	4.0	3.24	2.80
5	Δ-3-Carene	0.9992	0.32	0.98	3.7	3.85	-
6	α-Terpinene	0.9991	0.36	1.08	3.5	3.95	-
7	d-Limonene	0.9995	0.25	0.76	3.9	3.47	-
8	P-Cymene	0.9992	0.33	1.01	3.6	5.71	-
9	Eucalyptol	0.9999	0.11	0.33	3.8	4.54	-
10	Ocimene	0.9998	0.15	0.54	3.5	4.59	-
11	¥-Terpinene	0.9996	0.23	0.71	3.7	6.08	-
12	Terpinolene	0.9990	0.36	1.10	3.5	5.89	-
13	Linalool	0.9980	0.51	1.55	3.7	5.43	6.84
14	Isopulegol	0.9983	0.48	1.46	4.0	5.88	-
15	Geraniol	0.9964	1.21	3.66	4.3	4.67	-
16	β-Caryophllene	0.9992	0.33	1.00	4.0	7.02	-
17	α-Humulene	0.9985	0.46	1.38	3.5	6.38	4.53
18	Nerolidol-1	0.9969	1.31	3.96	3.6	3.13	-
19	Guaiol	0.9982	0.84	2.55	3.3	4.97	8.95
20	Caryphyllene oxide	0.9999	0.28	0.85	3.6	2.09	-
21	α-Bisabolol	0.9935	1.98	5.99	4.4	2.63	38.1

CONCLUSION

The Scion 4X6-GC analyser equipped with a split/spitless injector, Scion Instruments column, SQ-MS and HT3 Headspace Sampler is capable of analysing Terpenes from cannabis products in a qualitative and quantitative way.

The method developed is well suited for Terpenes analysis. It should be taken into account that the concentration of Terpenes in the sample can be very low and that a increase of sample volume or enrichment is necessary.

SCION Instruments

UK

Livingston Business Centre Kirkton South Road, Livingston West Lothian EH54 7FA Scotland, UK Phone + 44 1506 300 200 sales-eu@scioninstruments.com

The Netherlands

Amundsenweg 22-24 4462 GP Goes, The Netherlands Phone + 31 (0) 113 287 600 sales-eu@scioninstruments.com