



Comparison of BTEXS in Olive Oils by Static and Dynamic HT3

KEY WORDS: BTEX, Olive oil, MS, Flavour

INTRODUCTION

The health benefits of consumption of olive oils as part of a healthy diet reaches back to the mid 1950's. In the late 1990's, interest in characterizing other components of olive oils included determining the concentration of benzene, toluene, ethylbenzene, xylenes and styrene (BTEXS). Some of these compounds are naturally occurring, while others may be introduced through the processing of olive oil.

The method suggested for the determination of BTEXS in olive oil is static headspace. This application note will reproduce some of the static parameters found in this suggested method with the static features of the SCION Instruments HT3.

This application note will also demonstrate the versatility of the static and the dynamic capabilities of the HT3. These features allow laboratories the capability to use both static and dynamic headspace methods to characterize the volatile flavour components of olive oils and accurately quantify low level compounds like BTEXS.

Figure 1. SCION HT3 Headspace Sampler together with the SCION Instruments 8300 GC platform in combination with to 8700 SQ-MS equipped with 8400PRO autosampler .



Comparison of BTEXS in Olive Oils by Static and Dynamic HT3

INTRODUCTION

Olive oil has been part of the Mediterranean culture for thousands of years. Ancel Keys of the University of Minnesota was intrigued by the low incidence of coronary heart disease in the adult population of Naples Italy in 1952. His research has led in part to the current interest in the Mediterranean diet. This diet includes using olive oil, an oil low in saturated fats typically found in animal fats, palm, and coconut oil, and also a rich source of monounsaturated fats. Since then, other healthy uses for olive oil, such as salad dressing, have been discovered.

Because of the natural origin of olive oil, recent methods have been used to determine the amounts of benzene, toluene, ethylbenzene, xylenes and styrene (BTEXS) present. These methods recommend static headspace techniques to extract the BTEXS from various olive oils. These methods use a variety of techniques including differing sample sizes, GC/MS with selected ion monitoring (SIM), and stirring the sealed sample with Teflon coated stir bars to achieve the lowest detectable levels of these compounds.

This application note will explore the ability of the static and dynamic methods of the Teledyne Tekmar HT3 to detect BTEXS in various olive oils. The static method is utilized to characterize volatile flavor compounds responsible for the flavor and aroma of olive oils. The dynamic method is utilized for the concentration of low concentration of volatile organic compounds, like BTEXS, and to enhance the signal to noise (s/n) ratio for lower detection levels. The HT3 allows both methods to be performed in the same sequence.

Seven olive oil samples were obtained from a local market. Three of these were from the same manufacturer and consisted of a light flavour, pure, and extra virgin olive oil. The four remaining olive oils were all extra virgin olive oils from different manufacturers, including one labelled as organic. Six of the olive oils were packaged in glass containers. One was packaged in a plastic container. The labels of all of the olive oil indicated countries of origin.

Table 1. Instrumentation operating conditions.

Injector	Split 15 ml/min, 200 °C	
Column	SCION-624MS	
Oven Program	35°C (2 min), 10°C/min to 130°C, 20°C/min to 240°C (2 min)	
Carrier	Helium	
Flow	1.5 ml/min	
Autosampler	8400PRO	
Software	Compass CDS/ MSWS/ HT3 Teklink	
MS transfer line temperature	230°C	
Ion source temperature	230°C	
Ionization mode	EI	
Scan start	0.5	
Scan mode	Full Scan	
HT3	Static Value	Dynamic Value
Oven temperature	140°C	170°C
Transfer line temperature	140°C	180°C
Sample temperature	90°C	90°C
Sample equilibrium time	20 min	10 min

STANDARD SAMPLE PREPARATION

Seven different olive oil samples were obtained from a local market. Three of these, identified as S1, were from the same manufacturer and consisted of a light flavour, pure, and extra virgin olive oil.

The four remaining olive oils were all extra virgin olive oils from different manufacturers, including one labelled as organic. Six of the olive oils were packaged in glass containers. One, labelled S2, was packaged in a plastic container. The labels of all of the olive oil indicated countries of origin. Table 2 lists the label ingredients.

Comparison of BTEXS in Olive Oils by Static and Dynamic HT3

Table 2. Label statements of the seven locally available olive oils indicating the package container, olive oil grade, and the source country if the oil.

Supplier	Package	Grade	Olive Oil Origin	Package Origin
S1	Clear Glass	Lighter Flavor	Italy	Italy
S1	Clear Glass	Pure	Italy, Spain	Italy
S1	Clear Glass	Extra Virgin	Spain, Tunisia	Italy
S2	Clear PET Plastic	Extra Virgin	Spain, Tunisia	
S3	Green Glass	Extra Virgin	Italy	
S4	Clear Glass	Extra Virgin	Italy, Spain, Greece, Tunisia	
S5	Dark Green Glass	Organic Extra Virgin	Italy	

The samples for the static headspace analysis were prepared by placing 10g of each olive oil into separate 22mL headspace vials. The samples for the dynamic headspace analysis were prepared by placing 2g of each olive oil into separate 22mL headspace vials. All headspace vials were sealed with Teflon lined septa and capped and sealed with aluminium crimp caps.

An environmental standard containing benzene, toluene, ethylbenzene, o-, m-, and p-xylene and styrene was spiked into one of the olive oil samples to determine the retention times of these compounds. 1µL of a 2,000µg/mL solution was added to 10g of olive oil in a 22mL headspace vial.

All of the static and dynamic olive oil headspace vials were analyzed in a single sequence on the HT3 and the SCION Instruments GC/MS system with the parameters listed in Table 1.

Table 4. Normalized peak area for BTEXS observed in seven different olive oil samples by static and dynamic headspace methods.

Compound	Method	S1 Light	S1 Pure	S1 Ex Vir	S2	S3	S4	S5
Benzene	Static	100.0	35.8	13.5	14.4	18.3	12.8	15.7
	Dynamic	100.0	67.4	51.7	57.7	77.7	42.4	73.5
Toluene	Static	11.6	25.7	68.6	37.0	100.0	48.2	28.4
	Dynamic	6.9	22.5	65.8	42.2	100.0	48.0	32.0
Ethylbenzene	Static	0.0	34.0	36.5	23.0	100.0	50.1	42.1
	Dynamic	7.0	43.2	45.3	34.9	100.0	57.4	19.6
m-,p-Xylene	Static	1.8	21.0	46.4	32.5	100.0	55.7	14.1
	Dynamic	4.0	29.6	46.9	32.9	100.0	57.4	18.3
o-Xylene	Static	0.0	26.0	45.2	29.8	100.0	54.4	15.8
	Dynamic	1.7	31.0	41.2	29.1	100.0	52.0	20.3
Styrene	Static	0.0	4.4	100.0	17.1	15.3	51.8	10.7
	Dynamic	0.8	11.0	100.0	19.2	16.8	54.2	14.8

RESULTS

The Total Ion Chromatograms were used to quantitate the BTEXS in the samples using their primary quantitation ions typically used in environmental analysis. The primary quantitation ions and their secondary confirmation ions for the compounds are listed in Table 4.

Table 3. Primary quantitation ions and secondary confirmation ions for the BTEXS compounds.

Compound	Primary Quantitation ion m/z	Secondary Confirmation ion m/z
Benzene	78	77, 51
Toluene	91	92, 65
Ethylbenzene	91	106, 105
m-,p-Xylene	106	91, 105
o-Xylene	106	91, 105
Styrene	104	103, 78

The peak areas for the BTEXS were tabulated for the static and dynamic headspace methods. The peak areas were normalized to the largest peak area of the 7 olive oil samples for each headspace method. The normalized peak area data is presented in Table 4.

Comparison of BTEXS in Olive Oils by Static and Dynamic HT3

The BTEXS ions of 78m/z, 91m/z, 104m/z, and 106m/z were summed by the mass spectrometer software to create a summed ion chromatogram (SIC). The SIC of the static and dynamic headspace methods were compared for each olive oil sample. Figures 3 through 9 are these comparisons. The dynamic headspace method was at least 100 times more sensitive for these compounds in olive oil than the static headspace method.




Figure 2. Comparison of the Static and Dynamic Summed Ion Chromatograms of the BTEXS Ions for Olive Oil Sample S1 Lighter Flavor.




Figure 3. Comparison of the Static and Dynamic Summed Ion Chromatograms of the BTEXS Ions for Olive Oil Sample S1 Pure




Figure 4. Comparison of the Static and Dynamic Summed Ion Chromatograms of the BTEXS Ions for Olive Oil Sample S1 Extra Virgin




Figure 5. Comparison of the Static and Dynamic Summed Ion Chromatograms of the BTEXS Ions for Olive Oil Sample S2 Extra Virgin in PET Bottle




Figure 6. : Comparison of the Static and Dynamic Summed Ion Chromatograms of the BTEXS Ions for Olive Oil Sample S3 Extra Virgin.

Comparison of BTEXS in Olive Oils by Static and Dynamic HT3

CONCLUSION

Recent advances in characterization and enhanced detection limits have found trace amounts of BTEXS in samples of olive oils. One of the suggested methods for the determination of BTEXS in olive oils has been static headspace with GC/MS selection ion monitoring. The SCION Instruments HT3 is a flexible instrument that performs both static and dynamic headspace analysis for a wide variety of sample matrices including olive oils.

Seven olive oil samples were concurrently tested with both the static and the dynamic modes of the HT3 for benzene, toluene, ethylbenzene, xylenes and styrene. The static method requires 10g of olive oil compared to 2g for the dynamic mode.

The data indicated similar results for all of these compounds, except benzene. The dynamic mode detected higher concentrations of benzene in six of the seven olive oil sample. The seventh olive oil sample was used to normalize the peak area data.

The dynamic method also demonstrated increased signal to noise ratio for all of the compounds with 5 times less sample than the static method. The trap also allowed for the preliminary separation of the BTEXS compounds from other matrix effects from the olive oil.

The HT3 provides the olive oil industry the capability to perform both static and dynamic sampling of their products. This allows them to characterize olive oils for volatile flavour compounds in the static mode and then detect low levels of BTEXS in the dynamic mode.

ORDER INFORMATION

Part number	SCION HT3 Headspace
SC149300000	HT3 Headspace Autosampler 110 V. This static headspace autosampler comes with our 60 position autosampler, 10 position platen heater, integrated Optimix equilibrium system, and completely heated Silcosteel sample pathway. Entire system can be heated from ambient up to 300C in increments of 1C. Valve and loop configuration includes 1ml standard loop. System's gas flows and pressure are electronically monitored and controlled. Windows XP or 2000 is required to operate the included HT3 TekLink software. Power requirements are 115V, 50/60 Hz, 10A. Make sure to order the appropriate GC interface cable to ensure proper installation.
SC149300100	HT3 Headspace Autosampler 230 V. This static headspace autosampler comes with our 60 position autosampler, 10 position platen heater, integrated Optimix equilibrium system, and completely heated Silcosteel sample pathway. Entire system can be heated from ambient up to 300C in increments of 1C. Valve and loop configuration includes 1ml standard loop. System's gas flows and pressure are electronically monitored and controlled. Windows XP or 2000 is required to operate the included HT3 TekLink software. Power requirements are 230V, 50/60 Hz, 10A. Make sure to order the appropriate GC interface cable to ensure proper installation.
SC14930000S	HT3 Dynamic Headspace Autosampler 110V. This dynamic headspace autosampler comes with our 60 position autosampler, 10 position platen heater, integrated Optimix equilibrium system, and completely heated Silcosteel sample pathway. Entire system, including trap, can be heated from ambient up to 300C in increments of 1C. Valve and loop configuration includes 1ml standard loop. Dynamic mode includes 12" x 1/8" Tenax/Silica Gel and Charcoal (#3) trap and Vocab 3000 (K) trap. Static and Dynamic modes may be used interchangeably throughout the same schedule. System's gas flows and pressure are electronically monitored and controlled. Windows XP or 2000 is required to operate the included HT3 TekLink software. Power requirements are 115V, 50/60 Hz, 10A. Make sure to order the appropriate GC interface cable to ensure proper installation.
SC14930010S	HT3 Dynamic Headspace Autosampler 230V. This dynamic headspace autosampler comes with our 60 position autosampler, 10 position platen heater, integrated Optimix equilibrium system, and completely heated Silcosteel sample pathway. Entire system, including trap, can be heated from ambient up to 300C in increments of 1C. Valve and loop configuration includes 1ml standard loop. Dynamic mode includes 12" x 1/8" Tenax/Silica Gel and Charcoal (#3) trap and Vocab 3000 (K) trap. Static and Dynamic modes may be used interchangeably throughout the same schedule. System's gas flows and pressure are electronically monitored and controlled. Windows XP or 2000 is required to operate the included HT3 TekLink software. Power requirements are 230V, 50/60 Hz, 10A. Make sure to order the appropriate GC interface cable to ensure proper installation.

SCION Instruments

UK

4 Michaelson Square
Livingston
EH54 7DP, Scotland, UK
Phone +44 1506 300 200
sales-eu@scioninstruments.com

The Netherlands

Amundsenweg 22-24
4462 GP Goes,
The Netherlands
Phone +31 (0) 113 287 600
sales-eu@scioninstruments.com