

Plant Protection Products Impurity Screening by GC-FID with GC-MS Confirmation

AN165v1; SCION Instruments

Introduction

Impurities are an inevitable part of a chemical product; these are often process related, a by product of the synthesis or production of a product, or can be caused by degradation of the active ingredient (AI) or other components present within the product.

Determining impurity content is a requirement across many regulated industries from pharmaceutical, food, agriculture, and environmental.

But how do we decide which impurities are present within a product? Often chemical products can be very complex with dozens of impurities present, do we need to measure them all, or can we discount some, how do we decide this?

In this application note we will discuss how to determine which impurities should be identified and accurately quantified to satisfy regulatory requirements.

Regulations will vary by product type and by region, guidance documents will be available which will define the acceptable levels. SCION instruments recommends checking with local authorities in your region of interest before performing testing.

In this application note we will look at EU and UK regulations for submission of technical active substances and PPP No 283/2013 (Annex Section 4) and 284/2013 (Annex Section 5). These regulations have a guidance document SANCO/3030/99 rev5¹ which sets out the testing required for a successful submission of a product.

The product chosen for testing in this application was a commercially available bottle of technical grade active ingredient (TGAI) eugenol. Eugenol is a terpene which makes up around 80% of clove leaf oil and is attributed to the distinctive smell of cloves. Eugenol is also a powerful insecticide and when imported for this use falls under the requirements of the PPP regulations.²

The chemical structure and other details of eugenol can be found in Table 1.

There are two distinct types of impurity that must be considered. Relevant impurities which are defined as being "impurities of toxicological and/or ecotoxicological or environmental concern". These are defined by the product, for example eugenol, when submitted for use as a PPP, UK and EU regulations define methyl eugenol as a relevant impurity that can be present at a maximum of 0.1% of the technical material.³


The other type of impurity to consider is significant impurities. These are impurities that are present at or above a % defined by the relevant regulation. So for example for PPP products in the EU/UK the SANCO guidance document dictates that any impurity >0.1%.

Relevant impurities should be included in any impurity assay with a limit of quantification (LOQ) \leq the defined regulatory maximum, even if they are not found in the technical material.

Significant impurities are a little trickier, as these are often unknowns. Firstly we need to determine which impurities require identification and quantification. The simplest way to go about this is to use an analytical standard of the active ingredient – in this case eugenol. We can assume that any impurities present within the technical material are likely to be very similar in nature to eugenol. Therefore we can perform a screening of all present impurity peaks and compare this to a 0.1% Eugenol standard (impurity limit defined by SANCO).

Once we have these results we are then able to exclude any impurities that have been calculated to be <0.1%. The remaining impurities >0.1% will then need to be identified using mass spectrometry (MS) and finally their identity confirmed using the relevant analytical reference materials.

Table 1 Test Item Details

Eugenol TGAI	
IUPAC Name	4-Allyl-2-methoxyphenol
Molecular Weight	164.2 g/mol
Formula	C ₁₀ H ₁₂ O ₂
Label Purity Claim	>98%
Structure	

Impurity Screening

A nominal concentration of 10 mg/mL with regards to the technical material was chosen due to the suitable response shown by the impurities at this concentration.

A technical sample was prepared by weighing 100 mg eugenol TGAI (CoA claim 98% purity) into a 10 mL volumetric flask, the sample was made to volume with acetonitrile (HPLC grade, 98%) and mixed well.

Plant Protection Products Impurity Screening by GC-FID with GC-MS Confirmation



AN165v1; SCION Instruments

Acetonitrile blanks were ran to confirm that any impurities detected in the eugenol sample were not caused by impurities in the solvent.

A eugenol standard stock was prepared by weighing 50 mg eugenol analytical reference material (purity 99.8%) into a 100 mL volumetric flask, sample was made to volume with acetonitrile and mixed well. A dilution was then performed by pipetting 0.5 mL of this stock into a 25 mL volumetric flask, sample was made to volume with acetonitrile and mixed well to give 0.1% eugenol standard.

The technical sample was ran in duplicate, bracketed by the 0.1% standard, on an 8500 GC equipped with split/splitless (S/SL) injector and a flame ionization detector (FID). Method conditions can be found in table 2.

Table 2 GC-FID Instrument Parameters

Part	Settings
Autosampler	SCION 8400 PRO
Injector	300°C Split ratio 20:1
Injection Volume	1 µL
Column	SCION-5 30 m x 0.32 mm x 0.25 µm
Liner	SCION 1177 4 mm S/SL Focus
Carrier Gas	Nitrogen 1.5 mL/min
Oven Program	80°C (hold 1 min), 20°C/min to 300°C (hold 1 min)
Detector	Flame Ionization Detector (FID) Temperature: 325°C Air: 300 mL/min Hydrogen: 30 mL/min Make up (N ₂): 30 mL/min
Run Time	13 min
Software	CompassCDS

13 potential significant unknowns had been identified using this method. These unknowns were quantified using the average of the standard injections. Results can be see in Table 3.

The results show 3 impurities as being present ≥0.1%, unknowns 3, 7, and 11. These 3 unknowns will now need to be identified using MS.


An example chromatogram of the eugenol TGA1 sample can be seen in Figure 1.

Table 3 Screening Results for Eugenol TGA1

Impurity	% w/w
Unknown 1	0.02
Unknown 2	0.01
Unknown 3	0.15
Unknown 4	0.06
Unknown 5	0.05
Unknown 6	0.06
Unknown 7	2.70
Unknown 8	0.06
Unknown 9	0.04
Unknown 10	0.08
Unknown 11	0.12
Unknown 12	0.05
Unknown 13	0.05

Chromatogram: AN164_PPP_01_24_07_2023 13_36_59_Sample_2_1_channel1

System : AppGC
Method : AN165_PPP_Screening_01
User : User1
Acquired : 24/07/2023 13:40:40
Processed : 31/10/2023 12:14:30
Printed : 05/12/2023 15:47:08



Peak results :

Index	Name	Time [Min]	Area [µV.Min]	Tech
1	Unknown 1	4.35	11.5	
2	Unknown 2	4.40	2.3	
3	Unknown 3	5.46	77.0	
4	Eugenol	6.25	54619.9	
5	Unknown 4	6.29	30.2	
6	Unknown 5	6.49	24.7	
7	Unknown 6	6.72	28.9	
8	Unknown 7	6.88	1346.6	
9	Unknown 8	7.61	29.2	
10	Unknown 9	7.71	20.2	
11	Unknown 10	7.75	40.3	
12	Unknown 11	8.44	60.7	
13	Unknown 12	11.91	24.9	
14	Unknown 13	12.04	24.8	
Total			56341.3	

Figure 1 Example Chromatogram of a eugenol TGA1 sample

Plant Protection Products Impurity Screening by GC-FID with GC-MS Confirmation



AN165v1; SCION Instruments

MS Identification of Unknowns

In order to identify the unknown compounds a diluted eugenol TGA1 sample (1 mg/mL) was ran on an 8500 GC with 8700 MS single quad (SQ). Method conditions can be found in Table 4.

Table 4 GC-MS Instrument Parameters

Part	Settings
Autosampler	SCION 8400 PRO
Injector	240°C Split ratio initial 50:1 0.01 min off 0.30 min 50:1
Injection Volume	1 µL
Column	SCION-5MS 30 m x 0.25 mm x 0.25 µm
Liner	SCION 1177 4 mm S/SL Focus
Carrier Gas	Helium 1.0 mL/min
Oven Program	50°C (hold 1 min), 20°C/min to 300°C (hold 1 min)
Detector	Mass Spec (MS) Transfer Line 250 °C Scan 50 – 300 amu Solvent Delay 5 mins
Run Time	13 min
Software	MSWS with NIST spectral library

Using the powerful NIST spectral library search tool the MSWS software identified the most likely matches for each of the unknowns. The list of identified compounds can be seen in Table 5. For ease the 3 identified impurities shall be referred to as impurity 1, 2, and 3.

Example MSWS reports can be seen in Figures 2 and 3 for two of the unknown compounds.

Table 5 Identified Impurities

Unknown	Identity	CAS Number
3	4-(Prop-2-en-1-yl) phenol	501-92-8
7	3-(4-Isopropylphenyl)-2-methylpropionaldehyde	103-95-7
11	3-(4-hydroxy-3-methoxyphenyl)-2-propenal	458-36-6




Figure 2 MSWS Report for Unknown 3




Figure 3 MSWS Report for Unknown 7

Confirmation Using Analytical Standards

Analytical reference standards of the identified impurities were sourced in order to confirm the correct identity of each impurity. Once the identity of each impurity has been confirmed these standards can then be used in the validation and assay of an impurity method in order to support the successful submission of the product to the relevant regulatory bodies.

3-(4-hydroxy-3-methoxyphenyl)-2-propenal was unable to be sourced at the time of testing from a reliable and cost effective source. Therefore only the other 2 remaining impurities were sourced.

The EU defines the minimum purity required for Eugenol to be used as a PPP to be $\geq 99.0\%$ w/w. This particular product was imported under REACH regulations and could not be used as a PPP. We will therefore use the two impurities sourced as an example of analyte confirmation for a PPP.

In order to confirm the identities of the 2 impurities a eugenol TGA sample (prepared at 10 mg/mL in acetonitrile) was run under the conditions described in Table 2. Fortified eugenol TGA samples were also prepared (10 mg/mL in acetonitrile) and spiked with the relevant impurity (0.1% unknown 3 and 1.4% unknown 7). These samples were also run according to the method conditions in table 2.

If the impurities have been correctly identified we would expect to see an increase in the relevant peak areas of the fortified samples compared to the eugenol TGA sample.

Overlaid chromatograms of the eugenol TGA sample and the fortified eugenol TGA sample (spiked with impurities) are shown in figures 4 and 5.




Figure 4 Chromatogram showing overlay of peak for unknown 3 in a eugenol TGA sample and fortified eugenol TGA sample




Figure 5 Chromatogram showing overlay of peak for unknown 7 in a eugenol TGA sample and fortified eugenol TGA sample

The fortified samples show a clear increase in the peak area of the target analytes, with an exact match in retention time, this along with the data from the NIST spectral library search shows successful identification of the significant impurities present within the eugenol TGA batch.

The next stage after the significant impurities have been successfully identified would be a validation and assay of the relevant impurity (methyl eugenol) along with the identified significant impurities an example of this testing can be found in our application note AN166 Impurity validation and assay by GC-FID with hydrogen carrier gas.

A method validation and assay should also be performed for the active ingredient (AI), an example of this testing can be found in our application note AN164 plant protection products active ingredient assay GC-FID with nitrogen carrier gas.

Plant Protection Products Impurity Screening by GC-FID with GC-MS Confirmation

AN165v1; SCION Instruments

Conclusion

Impurities are an important part of regulatory testing across a variety of industries. In this application note we have seen how to screen for impurities in a plant protection product or pesticide in accordance with EU regulations.

Regulations will define relevant impurities for a product and define a limit above which each significant impurity must be identified and quantified in order for successful submission of the product for use under this particular regulation.

SCION instruments recommends checking with local regulatory authorities to ensure all testing and reporting requirements are met, or contact the SCION applications team for assistance.

The principles used in this application note can be applied to not only insecticides, but to other environmental pesticidal products such as biocides and fungicides as well as across other industries such as cosmetics, food, and drug testing.

Through quantification against an analytical reference standard of the active ingredient (eugenol) we were able to determine which impurities were significant. Then by using GC-MS in tandem with the powerful NIST library search tool we were able to identify these components. Finally the identity of the impurities was confirmed using analytical reference standards.

This application note is part of a series of 3, for active ingredient method validation and assay see application note AN164, for impurity method validation and assay see application note AN166.

For more information, please contact:

T(UK): +44 (0) 1506 300 200

T(EU): +31 (0) 113 287 600

E: sales-eu@scioninstruments.com

W: www.scioninstruments.com

Ordering Information

Ordering Information for the 8300 GC	
Part	Part Number
FID with Electrometer, DEFC Type 11, 120V	4362100101F
FID with Electrometer, DEFC Type 11, 230V	4362100102F
8400 PRO Autosampler for 8300 and 8500 GC	840000001
Suggested Consumables	
Part	Part Number
SCION-5 30 m x 0.32 mm x 0.25 µm GC Column	SC30233
SCION-5MS 30m x 0.25mm x 0.25µm	SC32223
Liner 1177 4MM S/SL FOCUS PK/5	41312101
15% Graphite/85% Vespel Ferrule 1/16" with 0.5 mm hole pk/10	41312149
BTO Septa 9 mm, pk/50	CR298713
10 µL fixed needle syringe, 5 cm, 0.47 mm OD, 26 g conical needle	41312133
Vial, 9-425 Screw Thread, 2 mL Clear Glass 12 x 32 mm Flat Base with Label, pk/100	41311000
Cap, Screw, Blue 9-425 Open Top Ribbed with 9mm Red PTFE/White Silicone Septa 1mm Thick. 100pcs/pk.	41311002

For ordering info on the SCION 8500 GC, which offers greater functionality with the option of up to 4 detectors (including MS), please contact your local SCION sales representative

References

1. European Commission SANCO/3030/99 rev.5, 22 March 2019 https://food.ec.europa.eu/system/files/2019-03/pesticides_ppp_app-proc_guide_phys-chem-ana_3030.pdf
2. Commission implementing regulation (EU) No 546/2013, Approval of Eugenol for use as a PPP, 14 June 2013 lex.europa.eu/LexUriServ/LexUriServ.do?uri=OJ:L:2013:163:0017:0020:EN:PDF
3. Commission Implementing Regulation (EU) NO 546/2013 14 June 2013 <https://www.legislation.gov.uk/eur/2013/546/2020-01-31?view=plain>