

Performance of EPA Method 8270 using Hydrogen Carrier Gas on a SCION GC-MS

AN0018 SCION Instruments

Introduction

United States Environmental Protection Agency (USEPA) Method 8270 is an analytical method for the detection of semi-volatile organic compounds in solid waste matrices, soils, air sampling media and water samples, by gas chromatography with mass spectrometry (GC-MS). The method measures a mixture of acids, bases, and neutrals in sample extracts. The complexity of these extracts demand a robust instrument that is easy to operate and maintain. Adding to method complexity is the uncertainty in both cost and supply of helium, forcing laboratories to consider hydrogen as a carrier gas. Hydrogen is not an inert gas; it is reactive and can be an explosion hazard if allowed to build up in the GC oven or manifold of the MS.

The SCION helium free analyser will ensure safe routine operation, with no performance change when operating under EPA Method 8270 specifications. SCION's unique axial ion source provides excellent robust operation and minimises unwanted protonation and spectral distortions. In addition, the GC with split/splitless (SSL) injector and inert pathway prevent compound degradation and reactions with the hydrogen carrier gas. This application note demonstrates the exceptional performance of the SCION GC-MS when operated under Method 8270 specifications.

Experimental

Table 1 details the operating parameters used throughout this application note. Please note, when using hydrogen carrier gas for the first time, more hydrocarbon background will be observed, but this will reduce with operation. The initial background can be reduced significantly by increasing the ion source temperature to 350°C, hydrogen column flow at 4mL/min and filament turned on for four hours.

Calibration standards containing 76 target compounds were used. The calibration ranged from 1ppm to 200ppm for the majority of compounds. Internal and surrogate standards were added at a concentration of 40ppm, in dichloromethane.

A pulsed split injection was used to minimise contact and residence time of compounds in the inlet. This is critical in hydrogen carrier gas due to its low viscosity and tendency to react with dichloromethane and form HCL.

The single goose-neck 4mm open inlet liner is preferred with Method 8270. As the liner does not contain glass wool, compound degradation is eliminated especially when using hydrogen as a carrier gas, due to its reactive properties.

Table 1. Analytical conditions of the GC-MS

Conditions	
S/SL	290°C, pulsed split 0.5µL (40psi for 0.3min, 70mL/min)
Liner	Single goose neck 4mm open inlet
Column	SCION-5MS 20m x 0.18mm x 0.18µm
Oven Programme	45°C (hold 3 min), 25°C/min to 100°C (1 min), 10°C/min to 310°C
Carrier Gas	Hydrogen 1mL/min
MS	SCION Select Helium Free
Scan Range	45-500Da
Ion Source	330°C
Transfer Line	300°C

Results

In order for hydrogen to be used as a carrier gas, the specifications of Method 8270 must be met. These specifications include tuning of the instrument, resolution, calibration, peak shape (Gaussian Factor) and system performance checks (SPCCs). Additionally, dichloromethane must be used as solvent to minimise degradation in the inlet. The GC-MS system must also produce mass spectra that match the NIST library and demonstrates robust operation when heavy matrices are analysed.

The SCION analyser can be auto-tuned with a tune-to-target feature for decafluorotriphenylphosphine (DFTPP), as required by the EPA method. Figure 1 shows the spectra for DFTPP acquired during the tune report. The set criteria for the DFTPP tune can be found in table 2. All requirements were met and the tune passed inspection. The DFTPP concentration was 50µg/mL with a 1µL injection.

Table 2. Tune acceptance criteria and true values

m/z	Acceptance Criteria	Value
51	30-60% of m/z 198	46.95
68	<2% of m/z 69	0.02
69	Present	40.10
70	<2% m/z 69	1.44
127	40-60% of m/z 198	40.96
197	<1% of m/z 198	0.33
198	base peak	100
199	5-9% of m/z 198	6.93
275	10-30% of m/z 198	21.40
265	>1% of m/z 198	2.34
441	<mass% of m/z 442	87.90
442	>40% of m/z 198	68.94
443	17-23% of m/z 442	17.02

Once the system has been tuned, a calibration series was analysed. EPA Method 8270 implements a calibration range of 0.5ppm to 160ppm, however, at the request of a laboratory that require EPA Method 8270 at a level from 0.075ppm to 30ppm, a calibration at this range was also analysed. These low concentrations show the excellent capability of the system to

Performance of EPA Method 8270 using Hydrogen Carrier Gas on a SCION GC-MS

AN0018 SCION Instruments

detect the compounds lower than those tested by EPA. As recommended by the method, internal standards were used at a concentration of 40ppm. The SCION SQ Mass Spectrometry software has a unique feature, Compound Based Scanning (CBS), in which SIM ions for compounds are stored in a library. The scan information, compound retention times and individual dwell times are all stored and are easily selected and loaded directly into a data acquisition method allowing both simultaneous mixed scan mode.




Table 3 details the retention times and SIM ions used throughout this application.

Table 3. Retention time and SIM Ions of method compounds

Compound Name	Retention Time	SIM Ion
N-Nitrosodimethylamine	3.639	74
2-Fluorophenol	5.177	112
Phenol-d5	6.106	99.1
Phenol	6.119	94.1
Bis(2-chloroethyl)ether	6.212	93
2-Chlorophenol	6.290	128
1,4-Dichlorobenzene-d4	6.542	152
N-nitroso-di-n-propylamine	7.071	70
Nitrobenzene-d5	7.312	82
2,4-Dimethylphenol	7.892	107
2,4-Dichlorophenol	8.237	162
Naphthalene	8.524	128
4-Chloro-3-methylphenol	9.597	107
2-Methylnaphthalene	9.867	142.2
1-Methylnaphthalene	10.071	142.2
2,6-Dimethylnaphthalene	10.248	156
2,4,6-Trichlorophenol	10.491	196
2,4,5-Trichlorophenol	10.563	196
2-Fluorobiphenyl	10.673	172
2,3,4-Trichlorophenol	10.698	196
Biphenyl	10.885	154
Dimethylphthalate	11.633	163

Acenaphthylene	11.841	152.2
Acenaphthene	12.242	153.2
2,3,5-Trimethylnaphthalene	13.139	170
Diethylphthalate	13.278	149.1
Fluorene	13.444	149.1
Azobenzene	13.853	77
2,4,6-Tribromophenol	14.023	330
Hexachlobenzene	14.731	284
Pentachlorophenol	15.240	266
Phenanthrene-d10	15.658	188.1
Penanthrene	15.714	178.2
Anthracene	15.834	178.2
Carbazole	16.281	167
1-Methylphenanthrene	17.188	192
Di-n-butylphthalate	17.195	149.1
Fluoroanthene	18.555	202.2
Pyrene	19.089	202.2
Terphenyl-d14	19.544	244.3
Butylbenzenepthalate	20.784	149.1
Chrysene	21.975	228.2
Benzo(a)anthracene	22.062	228.2
Bis(2-ethylhexyl)phthalate	22.176	149.1
Di-n-octylphthalate	23.649	149.1
Benzo(b)fluoranthene	24.376	252.2
Benzo(k)fluoranthene	24.436	252.2
Benzo(e)pyrene	24.945	252.2
Benzo(a)pyrene	25.047	252.2
Perylene-d12	25.171	264.2
Perylene	25.224	252.2
Indeno(123-cd)pyrene	27.469	278.3
Diben(ah)anthracene	27.531	278.3
Benzo(ghi)perylene	28.149	276.2

The calibration data, comprising of correlation coefficient, relative response factor and relative standard deviation, for calibration check compounds (CCC), as specified in Method 8270, can be found in Table 4. All quality control and system performance check samples passed the method requirements.

Table 4. Calibration data for SPCCs and CCCs

QC Compounds	Corr. Coeff.	RRF Value	RSD%
Phenol	0.9980	0.772	11.9
1,4-dichlorobenzene	0.9985	0.980	7.6
N-nitroso-di-propylamine	0.9991	0.719	8.5
2-nitro-phenol	0.9968	0.116	13.3
2,4-dichlorophenol	0.9981	0.114	8.4
hexachlorobutadiene	0.9991	0.315	7.0
4-chloro-3-methylphenol	0.9958	0.147	9.8
Hexachlorocyclopentadiene	0.9998	0.178	3.9
2,4,6-trichlorophenol	0.9997	0.142	8.0
Acenaphthene	0.9997	1.16	4.9
2,4-dinitrophenol	0.9993	0.070	14.5
4-nitrophenol	0.9940	0.134	6.8
pentachlorophenol	0.9994	0.082	14.2
Di-n-octyl phthalate	0.9998	0.199	4.6
Benzo(a)pyrene	0.9993	1.57	14.8

A biodiesel sample containing all 76 of the specified compounds were ran. The overall RSD% was 8.5%. The calibration curves of hexachlorocyclopentadiene, a system performance check

Performance of EPA Method 8270 using Hydrogen Carrier Gas on a SCION GC-MS

AN0018 SCION Instruments

compound, and pyridine, an EPA specified compound, can be seen in Figures 2 and 3.




Figure 2. Calibration curve of hexachlorocyclopentadiene




Figure 3. Calibration curve of pyridine

Figure 4 shows the calibration response factors for nine active compounds analysed. All response factors were less than 10% showing minimal reactivity when hydrogen is used as a carrier gas and dichloromethane as a solvent.




Figure 4. Calibration response factors of active compounds

When using hydrogen as a carrier gas, spectral quality can be compromised due to unwanted protonation or other reactions in the ion source. However, the Scion axial ion source, helium free package, minimises these reactions allowing good quality library matches to NIST. During this study, all seventy six compounds were detected using automated library search against NIST. Figure 5 illustrates an example library match for 1,2,4-trichlorobenzene with a match of over 96%.




Figure 5. Spectra match of 1,2,4-trichlorobenzene

EPA Method 8270 performance was tested using repeated injections of a contaminated sludge extract. A total of fifty injections were made, with the injection of the continuing calibration check (CCC) after every ten injections. The calculated concentration of the surrogates (target 40ppm) and the % recoveries in responses of the compounds were compared to the initial injection. Figure 6 shows the percent differences observed for the CCs.




Figure 6. Recovery % of calibration check compounds

APPLICATION NOTE

Performance of EPA Method 8270 using Hydrogen Carrier Gas on a SCION GC-MS



AN0018 SCION Instruments

The observed differences in concentrations throughout the fifty injections were considerably low. This demonstrates that the SCION SQ mass spectrometer is a robust instrument even when hydrogen is the carrier gas of choice.

Figure 7 shows the recovery of a surrogate standard over ten injections. The observed concentrations and % difference of a 40ppm 2,4,6-tribromophenol surrogate are also detailed.




Figure 7. Target and observed concentrations with % difference of 2,4,6-tribromophenol

Conclusion

EPA method 8270 is a challenging method covering a wide variety of compound classes and matrix types. As the cost and scarcity of helium rises, laboratories seek to amend methods to the use of hydrogen carrier gas. Due to its reactivity, hydrogen must be able to safely be used during GC-MS operation.

The SCION SQ has demonstrated excellent performance for method 8270, ensuring all qualitative and quantitative aspects of the method pass, and exceed, specification, when hydrogen is used as a carrier gas. The axial ion source and pulsed split injection technique, of the SSL injector, produced excellent library searchable mass spectra whilst also passing quality control criteria, even in heavy matrices. Additionally, the robustness of the ion source and sensitivity of the instrument allow lower reporting limits in challenging sample extracts, thus eliminating the challenges associated with EPA method 8270.

For more information, please contact:

T (UK): +44 (0) 1506 300 200

T (EU): +31 (0) 113 287 600

E: sales-eu@scioninstruments.com

W: www.scioninstruments.com