

SCION Sample Vials & Caps



SCION Sample Vials & Caps

2 ml Vials & Caps – Screw Top

- SCION Qualified vials & Septa for 8400/8410/CTC/EST samplers
- Narrow neck to ensure perfect septum closure
- White writing area for clear sample labeling
- Screw Cap allow for re-opening & closing vial
- Pre-slit septa to avoid septum damage

P/N	Description
41311000	Vial, 9-425 Screw Thread, 2mL Clear Glass 12x32mm Flat Base with Label. 100pcs/pk
41311001	Vial, 9-425 Screw Thread, 2mL Amber Glass 12x32mm Flat Base with Label. 100pcs/pk
41311002	Cap, Screw, Blue 9-425 Open Top Ribbed with 9mm Red PTFE/Sil Septa 1mm Thick. 100pcs/pk.
41311003	Cap, Screw, Blue 9-425 Open Top Ribbed with 9mm Red PTFE/Sil Septa 1mm Thick Pre-slit. 100pcs/pk.



SCION Sample Vials & Caps

2 ml Vials & Caps – Crimp Top

- SCION Qualified vials & Septa for 8400/8410/CTC/EST autosamplers
- Narrow neck to ensure septum closure
- White writing area for clear sample labeling
- Crimp Cap for perfect seals
- Pre-slit septa to avoid septum damage



P/N	Description
41311004	Vial, Crimp, Wide Opening, 2mL Clear Glass Flat Base 11mm with Label. 100pcs/pk.
41311005	Vial, Crimp, Wide Opening, 2mL Amber Glass Flat Base 11mm with Label. 100pcs/pk.
41311006	Cap, Crimp, Open, SILVER 11mm with Red PTFE/ White Silicone Septa 1mm Thick. 100pcs/pk.
41311007	Cap, Crimp, Open Top, CLEAR 11mm with Red PTFE/ White Silicone Septa 1mm Thick PRE-SLIT. 100pcs/pk.



SCION Sample Vials & Caps

10/20 ml Headspace Vials

- SCION Qualified vials & Septa for Tekmar/CTC/EST samplers. 10 ml fits 8410
- 10- or 20-ml size
- Magnetic/ non-magnetic caps
- Crimp or Screw top caps
- 10ml vials fit 8410 Autosampler



P/N	Description
41311008	Vial, Crimp, Headspace, 20mL Clear Glass 22.5x75mm. 20mm Beveled Edge. Rounded-Flat Bottom. 100pcs/pk.
41311009	Vial, Crimp, Headspace, 10mL Clear Glass 22.5x46mm. 20mm Beveled Edge. Rounded-Flat Bottom. 100pcs/pk.
41311010	Cap, Crimp, Open Top Aluminum, 20mm (10mm hole) with 20mm Natural PTFE/White Silicone Septa 3mm Thick. 100pcs/pk.
41311011	Cap, Crimp, Open Top Gold Magnetic, 20mm (10mm hole) with 20mm Natural PTFE/White Silicone Septa 3mm Thick. 100pcs/pk.
41311012	Vial, Screw Top Headspace, 20mL Clear Glass 22.5*75mm. 100pcs/pk.
41311013	Cap, Open Top Metal Screw, (8mm hole), 18mm Silver Color with 17.5mm Blue PTFE/White Silicone Septa 1.5mm Thick. 100pcs/pk.



SCION Sample Vials & Caps

40 ml EPA/TOC Vials

- SCION Qualified vials & Caps with septa for Purge & Trap
- 40 ml with 3mm septum for a 100% no worries perfect seal.
- Packaged includes vial & cap pre-assembled with dust cover to protect vials from dust and volatiles entering before sampling.
- 72 pcs per box, vials with caps.

P/N	Description
41311014	Cap, Screw, Open Top PP, 40mL 27.5x95mm Clear Glass EPA/TOC Vial 24-400 White with 22mm Natural PTFE/White Silicone 3.0mm thick Septa (EPA Quality) and Dust Cover. 72pcs/pk.



SCION Sample Vials & Caps

Crimp Tools

- The ideal crimp tools for 2-, 10- and 20-ml vials
- Cost effective manual crimpers
- Leak tight seal every time, for maximum sample integrity.

P/N	Description
99778700	Crimper for 11-mm crimp caps (2-mL vials)
99780400	Decapper for 11-mm crimp caps (2-mL vials)
392612019	Crimper for 20-mm crimp caps (10/20-mL vials)
190003617	Decapper for 20-mm crimp caps (10/20-mL vials)
