


# SCION Cold On Column (1093) Parts



#	Description	P/N
1a	Injector Nut	394966601
	Injector Nut Wrench	390842300
2	Automatic Start Switch	390820601
3	Septum, 7/16"	See Septa
	BTO, lowest bleed	CR298777
	Marathon, autosampler	CR239287
	Advanced Green 3	CR246725
4	Septum Pick	7200008400
	Septum Support, gold plated	391821100
5	Liner Insert	See Liners
	Glass Insert, default	SG092034
	High performance Insert	190010906
	On-column Insert	190010907
	Packed	392611950
6	Wave Spring Washer	1492000500
7	Coolant Exit Options	
	Frit Restrictor	3600004600
8	Plug, 1/8"tube	2824707100
	Screw	391866306
9	Injector Body	
10	Column Ferrule	See ferrules
	40%G/60%V, up to 400°C	CR213400
	Graphite, up to 450°C	CR211400
11	Column Nut	
	0.9mm Capillary, brass	394955100
	1.6mm for 1/16", brass	CP742351
	High Temp, stainless steel	CP743117
	Carrier Gas Pencil Filter, duo-bed	CP742210