


SCION FID Spare Parts



#	Description	P/N
1	Signal Cable	391783000
2	Signal Probe	200187300
3	Ignitor Cable	392548701
4	Ignitor Probe	200187200
5	Collector Tube	394958700
6	Lower Ceramic Insulator	2100003200
7	Screws,8-32x1/4"	1214200804
8	Tower Screws,8-32x3"	391866302
9	Tower Top Plate	3400073100
0	Detector Tower Body	100099300
11	FlameTip	
	0.01"capillarytipw/nut	200187500
	0.02"standardtipw/nut	200193800
12	FlameTip Ferrule	
	Vespel 100%,upto350°C	CR212200
	40%G/60%V,upto400°C	CR213200
	Graphite,upto450°C	CR211200
13	Aluminium Crunch Washer	1500334701
14	Detector Base	392547501
15	Column Ferrule	See Ferrules
16	Column Nut	
	0.9mm Capillary, brass	394955100
	1.6mm for1/16", brass	CP742351
	High Temp, stainless-steel	CP743117
	Probe O-Ring,0.505OD	390596004
	Flow Tube,1/8"OD 200196700	200196700
	FID TestMix	8200504807