


SCION Flash (1061) Spare Parts



#	Description	P/N
1a	Injector Nut, Autosampler	390812700
1b	Injector Nut, manual	392595501
	Injector Nut Wrench	390842300
2	Septum, 9.5 mm	See Septa
	BTO, lowest bleed	CR298705
	Marathon, autosampler	CR239188
	Advanced Green 3	CR246124
	Septum Pick	7200008400
3	Automatic Start Switch	390820601
4	Spring	Call
5	Cover & insulation	Call
6	Glass Insert	392611943
	Flash Glass Insert Megabore 5pk	392611943
	Flash Glass Insert Packed 5pk	392611944
	Packed Colum Adapter Kit 1/8" SS	392558892
7	Injector Body	392548301
8	Column Ferrule	See ferrules
	40%G/60%V, up to 400°C	CR213400
	Graphite, up to 450°C	CR211400
9	Injector Body Nut, stainless steel	SWSS4021
10	Flash Vaporization Column Guide	392558301
11	Column Ferrule	See ferrules
12	Column Nut	
	0.9mm Capillary, brass	394955100
	1.6mm for 1/16", brass	CP742351
	High Temp, stainless-steel	CP743117
	Carrier Gas Pencil Filter, duo-bed	CP742210