


# SCION PFPD Spare Parts



#	Description	P/N
1	Hex screw,8-32x3/4"	392519401
2	Ignitor Cap Assy	392513400
3	Ignitor Cable Assy	392546101
4	Screw,8-32x11/2"	391866307
5	Aluminium Crunch Washer	1500334701
6	Quartz Combustor	
	2 mm, Sulphur	392517600
	3 mm, Phosphorous	392517700
7	PFPD Tower	
	N,P,S Assy	392513301
	S and P Assy	392513300
8	Combustor Holder	
	2 mm, Sulphur	392517800
	3mm, Phosphorous	392517901
9	Aluminium Crunch Washer	1500334701
10	Seal Combustor Support	392513800
11	Detector Base	392547601
12	High Voltage Cable Assy	
13	Signal Cable Assy	
14	Column Ferrule	See Ferrules
15	Column Nut	
	0.9mm Capillary, brass	394955100
	1.6mm for1/16", brass	CP742351
	High Temp, stainless-steel	CP743117
	Diffusion Source kit	392515700
	Combustor Support Wrench	392519200
	PFPD TestMix	8200504803