


SCION PTV Inlet (1079) Spare Parts



#	Description	P/N
1a	Injector Nut	392595401
1b	Injector Nut	394966601
	Injector Nut Wrench	390842300
2	Automatic Start Switch	390820601
3	Septum, 11.5 mm	See Septa
	BTO, lowest bleed	CR298777
	Marathon, autosampler	CR239787
	Advanced Green 3	CR246725
	Septum Pick	7200008400
4	Septum Support, gold plated	391867600
5	Liner Ferrule, graphite 5mm	392534201
6	Liner	See Liners
7	Injector Body	
	EFC21, stainless-steel	392544001
	EFC21, UltiMetal™	392544011
	EFC25, stainless steel, manual	392559601
	EFC25, UltiMetal™, manual	392597304
8	Column Ferrule	See ferrules
9	Column Nut	
	0.9mm Capillary, brass	394955100
	1.6mm for 1/16", brass	CP742351
	High Temp, stainless steel	CP743117
	Split Vent Pencil Filter, charcoal	CP742211
	Carrier Gas Pencil Filter, duo-bed	CP742210