

# SCION Septa

---



# SCION Septa

## SCION BTO Septa

- Bleed & temperature optimized BTO Septa; the choice for use with Mass Spec
- Ultra low Bleed; Maximum operational temperature of 450°C.
- Centerguide for easy penetration and minimized septum damage & artifacts



P/N	Description	Diam. (mm)	Center Guide
CR298713	SCION BTO Septa with CenterGuide, 9 mm, pk/50.	9	Yes
CR298705	SCION BTO Septa 3/8", pk/50	9.5	No
CR298777	SCION BTO septa with centerguide 11.5mm, pk/50	11.5	Yes

## SCION Marathon Septa

- Very low bleed, softer material, gentle on autosampler needles
- Longer lifetime; up to 400 injections without fail with autosampler.
- Maximum operational temperature of 350°C.



P/N	Description	Diam. (mm)	Center Guide
CR239778	SCION Marathon Septa with CenterGuide, 9 mm, pk/50.	9	Yes
CR239188	SCION Marathon Septa 3/8", pk/50	9.5	No
CR239287	SCION marathon Septa with centerguide 11.5mm, pk/50	11.5	Yes

## SCION Advanced Green Septa

- Considerably longer lifetime, reduced injection port sticking, and low bleed.
- Softer material, gentle on autosampler needles
- Maximum operational temperature of 350°C.



P/N	Description	Diam. (mm)	Center Guide
CR246713	SCION Marathon Septa with CenterGuide, 9 mm, pk/50.	9	Yes
CR246124	SCION Marathon Septa 3/8", pk/50	9.5	No
CR246724	SCION marathon Septa with centerguide 11.5mm, pk/50	11.5	Yes