

LC6000 Promotion

LC Bundles with CompassCDS

The SCION 6000 series HPLC delivers confidence in results through outstanding lifetime performance. A robust design maximises uptime and productivity levels whilst minimising cost of operation.

With an array of automation options for workflow optimisations, the LC6000 is the chromatographer's choice in HPLC. The 6000 Series HPLC line-up comes with a portfolio of detectors including UV, DAD, FL and RID, all designed for stability, precision and stability.

All LC bundles come with CompassCDS, an industry-proven powerful networked chromatography data system software. With its intuitive, easy to use, information rich user interface CompassCDS provides operators with a powerful tool for a variety of operations in different analytical fields. Designed for complete instrument control, data acquisition, processing and reporting, CompassCDS is fully customisable providing a powerful data analysis tool with comprehensive calibration and calculation options, whilst being 21 CFR 11 compliant.




21 CFR 11 COMPLIANT

21 CFR 11 Provides full audit trail for superior data integrity.



100% SCALABLE

from Standalone to client server C/S solution



SCION Instruments

Livingston Business Centre
Kirkton South Road, Livingston
West Lothian EH54 7FA
Scotland, UK

Phone +44 (0)1506 300 200
sales-eu@scioninstruments.com

LC6000 Promotion

LC Bundles with CompassCDS

1

Complete LC System with UV
Valid until 31st December 2021

6510 Organiser:	Extra-large solvent cabinet 6 x 1L mobile phase bottles 3 x 500mL mobile phase bottles
6100 Pump:	Up to 600 bar maximum operation pressure Up to 10mL flow 6-channel degassing unit 4-channel low pressure gradient system
6220 Autosampler:	Peltier cooling system 4°C to -5°C ambient Direct injection system 200 x 1.5mL vials sample tray Microwell Plate optional
6310 Column Oven:	Peltier heating / cooling system Temperature range 1°C to 85°C Houses 3 x 250mm columns
6410 UV Detector:	Double beam photometric system 14µL flow cell Wavelength range 190nm to 600nm
Software:	CompassCDS
PC Hardware:	Precision Dell PC and Monitor
Extras Included:	C18 column (250mm x 4.6mm x 5µm)
Options:	IQ/OQ/ PQ* *Purchased separately

2

Complete LC System with DAD
Valid until 31st December 2021

6510 Organiser:	Extra-large solvent cabinet 6 x 1L mobile phase bottles 3 x 500mL mobile phase bottles
6100 Pump:	Up to 600 bar maximum operation pressure Up to 10mL flow 6-channel degassing unit 4-channel low pressure gradient system
6220 Autosampler:	Peltier cooling system 4°C to -5°C ambient Direct injection system 200 x 1.5mL vials sample tray Microwell Plate optional
6310 Column Oven:	Peltier heating / cooling system Temperature range 1°C to 85°C Houses 3 x 250mm columns
6430 DAD Detector:	13µL Flow cell Wavelength range 190nm to 900nm
Software:	CompassCDS
PC Hardware:	Precision Dell PC and Monitor
Extras Included:	C18 column (250mm x 4.6mm x 5µm)
Options:	IQ/OQ/PQ* *Purchased separately

3

Complete LC System with FL
Valid until 31st December 2021

6510 Organiser:	Extra-large solvent cabinet 6 x 1L mobile phase bottles 3 x 500mL mobile phase bottles
6100 Pump:	Up to 600 bar maximum operation pressure Up to 10mL flow 6-channel degassing unit 4-channel low pressure gradient system
6220 Autosampler:	Peltier cooling system 4°C to -5°C ambient Direct injection system 200 x 1.5mL vials sample tray Microwell Plate optional
6310 Column Oven:	Peltier heating / cooling system Temperature range 1°C to 85°C Houses 3 x 250mm columns
6440 FL Detector:	Flow Cell 12µL Ex Wavelength 200nm to 850nm Em Wavelength 250nm to 900nm
Software:	CompassCDS
PC Hardware:	Precision Dell PC and Monitor
Extras Included:	C18 column (250mm x 4.6mm x 5µm)
Options:	IQ/OQ/PQ* *Purchased separately

Contact us

Livingston Business Centre,
Livingston, EH54 7FA, Scotland, UK

Tel: +44 (0)1506 300 200
sales-eu@scioninstruments.com