


# Chromatogram: AN164\_PPP\_01\_17\_07\_2023 13\_17\_06\_14\_1\_channel1

System : APP\_GC\_04  
 Method : AN164\_PPP\_01  
 User : User1

Acquired : 17/07/2023 13:20:41  
 Processed : 19/07/2023 14:47:18  
 Printed : 19/07/2023 14:49:57



## Peak results :

Index	Name	Time [Min]	Quantity [% Area]	Height [µV]	Area [µV.Min]	Area % [%]
1	Eugenol	6.22	41.52	137894.6	2694.8	41.519
2	Internal Standard	7.63	58.48	179476.7	3795.8	58.481
Total						
			100.00	317371.3	6490.6	100.000

MP1